



SUPER SPORT™

Advantrode® Elite Electrode

- Specifically engineered ez-connect snap allows the *Super Sport*™ to be compatible with all Compex® brand units. (Edge, Performance, Elite and Wireless USA)
- Automated manufacturing process helps lead to a significant cost saving compared to Compex® brand electrodes.
- Layer of highly conductive 99.999% pure silver (Ag) on carbon film is proven to provide efficient dispersion of electrical current over the entire surface area of the electrode.
- Extra thick hydrogel adhesive provides maximum reusability throughout the life of the electrode.
- Strip of ultra-conductive film resting on bottom of snap helps further improve conductivity, while lowering impedance for optimum performance.
- Medical-grade blue spunlace for ultimate flexibility and conformability to body contours during therapy.



Item No: SSB22
2" X 2" sq. w/snap 4/pk



Made in the USA

PEPIN MANUFACTURING, INC.
1875 HWY 61 South, Lake City, MN 55041
800-291-6505 • 651-345-5655 • Fax: 651-345-5656
www.pepinmfg.com • Email: electrodes@pepinmfg.com

