

ADVANTRODE® ELITE

EXPERIENCE THE UNSURPASSED PERFORMANCE!

ADVANTRODE® ELITE electrodes are offered in a variety of flexible, durable topcoat materials: White Foam, Tan Tricot, Tan Spunlace, White Spunlace and Cool Blue. PMI's state-of-the-art manufacturing processes, superior medical-grade materials and commitment to quality is reflected in each electrode. The expanded selections in different sizes, shapes and hydrogels provide you and your customers the freedom to choose the most effective electrodes for maximum treatment results.

FEATURE

Highly conductive silver/carbon film

35 mil dual layer hydrogel

Tamper-evident bag
*(Available with most item numbers)

Copper stranded wire with tab and molded grip

Molded grip wire socket

Large selection of shapes and sizes

Foam and cloth topcoats

BENEFIT

Proven efficient dispersion of electrical current over the entire surface of the electrode.

Maximum reusability throughout the life of the electrode.

Assurance that the product has not been altered.

Enhanced current dispersion, lowest impedance. Tab prevents wire pull out.


Ease of connection and disconnection from lead wire.

Adaptable to multiple electrode application sites.

Ultimate flexibility and conformability to body contours.



AdvanTrode® Elite		Tan Spunlace	
Item #	Description	Elect/Pack	Pk/Case
WTS2	2" (5cm) Round	4	20
WTS3	3" (7.5cm) Round	4	20
WTS15	1.5" x 1.5" (4cm x 4cm) Square	4	20
WTS22	2" x 2" (5cm x 5cm) Square	4	20
WTS125o	1.5" x 2.5" (4cm x 6.5cm) Oval	4	20
WTS24o	2" x 4" (5cm x 10cm) Oval	4	20
WTS152	1.5" x 2" (4cm x 5cm) Rect.	4	20
WTS13	1.75" x 3.75" (4.5cm x 9.5cm) Rect.	4	20



AdvanTrode® Elite		White Foam	
Item #	Description	Elect/Pack	Pk/Case
WF1	1" (2.5cm) Round	4	20
WF375	1.375" (3.5cm) Round	4	20
SF375	1.375" (3.5cm) Round w/Snap	4	20
WF2	2" (5cm) Round	4	20
WF3	3" (7.5cm) Round	4	20
WF22	2" x 2" (5cm x 5cm) Square	4	20
WF125o	1.5" x 2.5" (4cm x 6.5cm) Oval	4	20
WF24o	2" x 4" (5cm x 10cm) Oval	4	20
WF13	1.75" x 3.75" (4.5cm x 9.5cm) Rect.	4	20
WF35	3" x 5" (7.5cm x 13cm) Rect.	2	20



AdvanTrode® Elite		White Spunlace	
Item #	Description	Elect/Pack	Pk/Case
WWS2	2" (5cm) Round	4	20
WWS3	3" (7.5cm) Round	4	20
WWS22	2" x 2" (5cm x 5cm) Square	4	20
WWS125o	1.5" x 2.5" (4cm x 6.5cm) Oval	4	20
WWS24o	2" x 4" (5cm x 10cm) Oval	4	20
WWS13	1.75" x 3.75" (4.5cm x 9.5cm) Rect.	4	20
WWS35	3" x 5" (7.5cm x 13cm) Rect.	2	20



AdvanTrode® Elite		Cool Blue	
Item #	Description	Elect/Pack	Pk/Case
EW15	1.5" x 1.5" (4cm x 4cm) Square	4	20
EW22	2" x 2" (5cm x 5cm) Square	4	20
EW13	1.75" x 3.75" (4.5cm x 9.5cm) Rect.	4	20

ELECTRODE SELECTION

SHAPES

ROUND

SQUARE

OVAL

RECTANGLE

BUTTERFLY

BACKLAP



AdvanTrode® Elite		Tan Tricot	
Item #	Description	Elect/Pack	Pk/Case
WT1	1" (2.5cm) Round	4	20
WT375	1.375" (3.5cm) Round	4	20
WT2	2" (5cm) Round	4	20
WT3	3" (7.5cm) Round	4	20
WT11	1.75" x 1.75" (4.5cm x 4.5cm) Square	4	20
WT22	2" x 2" (5cm x 5cm) Square	4	20
WT24o	2" x 4" (5cm x 10cm) Oval	4	20
WT35o	3" x 5" (7.5cm x 13cm) Oval	2	20
WT152	1.5" x 2" (4cm x 5cm) Rect.	4	20
WT13	1.75" x 3.75" (4.5cm x 9.5cm) Rect.	4	20
WT35	3" x 5" (7.5cm x 13cm) Rect.	2	20
WT47	4.25" x 7" (11cm x 18cm)	1	20
ST11	1.75" x 1.75" (4.5cm x 4.5cm) Sq. w/Snap	4	20
ST13	1.75" x 3.75" (4.5cm x 9.5cm) Rect. w/Snap	4	20



AdvanTrode® Elite		Back Therapy	
Item #	Description	Elect/Pack	Pk/Case
WF46	4" x 6" (10cm x 15cm) Butterfly	1	20
WF47	4.25" x 7" (11cm x 18cm) Rect.	1	20
WF58	5" x 8" (13cm x 20cm) Rect.	1	20
WF1315	13" x 1.5" (33cm x 4cm) Backlap	1	20
GVK	4.25" x 7" (11cm x 18cm) Rect.	1/Kit	10
(gohanic kit)	1.5" x 2.5" (4cm x 6.5cm) Oval	2/Kit	



AdvanTrode® Elite		Skin Friendly Blue Gel	
Item #	Description	Elect/Pack	Pk/Case
BWT11	1.75" x 1.75" (4.5cm x 4.5cm) Square	4	20
BWT13	1.75" x 3.75" (4.5cm x 9.5cm) Rect.	4	20
BWT35	3" x 5" (7.5cm x 13cm) Rect.	2	20
BST11	1.75" x 1.75" (4.5cm x 4.5cm) Sq. w/Snap	4	20
BST13	1.75" x 3.75" (4.5cm x 9.5cm) Rect. w/Snap	4	20